



accucold™
All-Freezer

PER OPTIMUM STORAGE CONDITIONS
Always keep the unit at least 20% full to maintain a stable temperature environment. If storing meat, fish, and frozen goods.

NOTE:
The evaporator fan stops on the fan sensor detecting an air temperature change. The temperature probe is located in the right drawer. All temperatures are monitored by the CCU.

ALARM/THERMOMETER
The alarm displays the temperature change in the drawers using a digital display. The temperature probe is located in the right drawer. All temperatures are monitored by the CCU.

DIGITAL THERMOSTAT
The digital thermostat allows you to adjust the temperature. The digital thermostat will automatically adjust the set point when the unit is in a temperature change.

Learn more at accucold.com



**Built-in undercounter
medical/scientific -25°C
capable all-freezer with
front control panel
equipped with a digital
thermostat and NIST
calibrated
thermometer/alarm**

SKU: VT65MLBIMED2

Price: \$2705

Category: [Lab Equipment & Supplies](#), [Refrigerator Freezers](#)